

CamAPS FX is a sophisticated insulin delivery app for people with type 1 diabetes that is very easy to use.

- CamAPS FX is an Android smartphone app that continually adjusts the insulin given via an insulin pump based on the glucose sensor data it receives.
 - It has to be connected via Bluetooth to the insulin pump and sensor to work, so the phone with the app needs to stay within 6m of the person.
- People using CamAPS FX will still need to carbohydrate count and give an insulin bolus, via the CamAPS FX bolus calculator, for any meals or snacks (ideally about 10-15 minutes before eating).

CamDiab Customer Support:

- Telephone support (call: 020 3695 3780) is available Monday-Friday 9-5pm
- Email support (support@camdiab.com)
- CamAPS FX user manual can be found on CamDiab's website (<u>www.camdiab.com/support</u>) or in the app's 'help' menu.
- Click on image to open user manual

CamAPS FX Online Training and Webinars:

- ✓ FREE online training and webinars are available here: www.camdiabtraining.com.
- ✓ Training available for
 - People with type 1 diabetes and their families
 - Diabetes healthcare professionals
 - Teachers and education support staff

CamAPS FX Data Monitoring for Parents and Carers:

- Create a personal Diasend account: <u>www.diasend.com</u>
- \checkmark Link the Diasend account to CamAPS FX main menu > share > user 1 Diasend
- Download the Diasend app on parent's or carer's phone
- Click on image for more info about the Diasend app.

Peer Support available via 2 closed Facebook groups:

- CamAPSFX for kids (for families using CamAPS FX)
- CamAPS users (for all CamAPS FX users)
- Click on the image to go to FB group

* Please note these are not CamDiab social media pages

CamAPS FX Social Media:

- Follow CamAPS FX on social media to keep in the loop...
- Click on the icons to follow.

Contact CamDiab:

Telephone: 020 3695 3780 Email: support@camdiab.com Website: www.camdiab.com













amDiab



